

# Christchurch Daily Equipment Check List

Christchurch Version 1

	DAILY CHECKS	Monday	Tuesday	Wednesday	Thursday	Friday	Saturday	Sunday
1	<b>Drug Checks</b> Two people (1 must be a clinician), to check these drugs <b>every day</b> . Midazolam, Diazepam & Morphine (Oramorph), Codeine NOT AVAILBLE Oxycodone, Morphine Sulphate NOT AVAILABLE	Evening Host	Evening Host	Evening Host	Evening Host	Evening Host	<b>Morning Host</b>	<b>Morning Host</b>
2	<b>Consulting Room Boxes (should be 5 including isolation room. )</b> Check contents against list on box, report to SM any missing items.1 2 5 6 7 8	Evening Host	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts
3	<b>Panic Alarms</b> All in a small box in host cupboard, ensure they are back in the box end of shift (6)	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts
4	<b>Patients and downloads folder (on desktop)</b> Please make sure its clear start and end of shift	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts
5	<b>Check respiratory hood and HCIS cleaning box are present and complete</b> Check contents against list on both boxes Clean contents and box itself with Clinell wipe	Evening Host	Evening Host	Evening Host	Evening Host	Evening Host	Evening Host	Evening Host
6	<b>Resus bag</b> - Location: Store cupboard in corridor Check contents as per contents sheet, replace any stock if needed. Check any OOD and log any imminent. Ensure green tick is showing which indicates Defib is working. Clean defib with Clinell wipe	Evening Host	Evening Host	Evening Host	Evening Host	Evening Host	<b>Morning Host</b>	<b>Morning Host</b>
7	<b>Blood Taking Box</b> Ensure missing items are replaced from stock, check for expiry dates	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts
8	<b>Sharp boxes</b> Should not exceed the 3-month date, should be closed but not locked, any overfilled/not closed boxes to be reported	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts	All Hosts
9	<b>Manual Scripts</b> Count and record manual prescription packs in Audit book	Overnight Host	Overnight Host	Overnight Host	Overnight Host	Overnight Host	Overnight Host	Overnight Host
10	<b>Emergency drugs</b> - Location: Inside drugs cupboard in store room Ensure drugs are in the appropriate drugs cupboard in store room	Evening Hosts	Evening Host	Evening Host	Evening Host	Evening Host	PM Host	PM Host
11	<b>Oxygen</b> - Location: Small and Large bottles Please accurately indicate in the signature box the levels as below: EMPTY / ¼ FULL / ½ FULL / ¾ FULL / FULL If empty, report to shift manager	Evening Host	Evening Host	Evening Host	Evening Host	Evening Host	Morning Host	Morning Host
	<b>Small cylindres</b>							
	<b>Large Cylinder</b>							
12	<b>Monitoring Drug Temperatures</b> – Location: Drug Cabinet Check temp and record Min and Max temp on form • If the alarms go off whilst on shift, please let the shift manager know	Overnight Host	Overnight Host	Overnight Host	Overnight Host	Overnight Host	Overnight Host	Overnight Host
		MIN: MAX:	MIN: MAX:	MIN: MAX:	MIN: MAX:	MIN: MAX:	MIN: MAX:	MIN: MAX:

# Christchurch Daily Equipment Check List

Christchurch Version 1

13	Handwash Audits	Eve	o/n	Eve	O/N	Eve	o/n	Eve	o/n	Eve	o/n	Am	Pm	o/n	Am	Pm	O/n	
	Log into Radar, check which clinicians need handwash. Complete survey. Initial when completed on Radar																	
	<b>Weekly Checks</b>	Monday		Tuesday		Wednesday		Thursday		Friday		Saturday		Sunday				
13	<b>Blood Taking Box</b>											Overnight						
	Check all expiry dates for stock and replace as required																	
14	<b>Health and Safety Checks Weekly</b>			Overnight														
	Complete H&S checklist, keep original in folder and flag any concerns																	
15	<b>Paediatric/Adult pulse oximeter - Location:</b> Store Room within Sats monitor box Weekly	Overnight																
	Switch on machine, place on finger and ensure it takes a reading. Ensure Adult, child and Infant leads are available Clean the machine and box itself with Clinell wipe																	
16	<b>Nebuliser Machine - Location:</b> Store room within the nebuliser box Weekly					Overnight												
	Switch on – does it make a noise? Clean the machine with Clinell wipe																	
17	<b>Doppler Weekly</b>											Overnight						
	Check Doppler is stored in cupboard and the equipment turns on. Clean the machine with Clinell wipe																	
18	<b>Blood Glucose Machine / Dual Keytone Weekly Check</b>							Overnight										
	Host to check machine as per guidelines Clean the machine and box itself with Clinell wipe																	